

## **SAMPLE:**

For XRD, Please provide fine powdered sample (around 25 mg for preferably) in a small vial.

For IR spectrometer, Please provide fine powdered sample (around 10 mg for preferably) in a small vial.

For UV spectrometer, Please provide fine powdered sample (around 10 mg for preferably) in a small vial.

For CV, Charge/discharge, Life Cycle, Impedance, EQCM and RRDE/RDA measurements, Please provide fine powdered sample (around 100 mg for preferably) in a small vial.

For Ball mill and Q Sonica for Nanoparticles, Please provide fine powdered sample (around 10 mg-5g for preferably) in a vial.